



Universal Tray



Summary

The DDO Universal Trays can be attached virtually to any wheelchair model. The trays provide support and encourage daily activities. There are two models of DDO Universal Trays:

Universal Wheelchair Tray – Hard Mounting

The Universal Wheelchair Tray – Hard Mounting comes in two sizes (Small for 12" and 14" wheelchair frame widths / Large for 16", 18" and 20" wheelchair frame widths)

The tray uses the wheelchair backrest vertical rails for mounting.

Universal Wheelchair Tray – Soft Mounting

The Universal Wheelchair Tray – Soft Mounting comes in two sizes (Small for 12" and 14" wheelchair frame widths / Large for 16", 18" and 20" wheelchair frame widths)

The tray uses the armrest pads for mounting.

Size

Size

Small

Suitable for 12" (30cm) and 14" (35cm) frame widths

Large

Suitable for 16" (40cm), 18" (45cm) and 20" (50cm) frame widths



Please contact us at info@diversable.org for more information

Please contact us at productfeedback@diversable.org to share your experience with our products in the field.

Lack of/or poorly fitted products can result in permanent deformity or health complications which can lead to death. A wheelchair provider with limited knowledge and skill can improve through direct training and supervised practice. Please note: The DDO Universal Tray Soft and Hard Mounting designs are the property of Diversability Development Organization. No part of these design, including all original concepts and ideas contained herein, may be reproduced or transmitted in any form, or by and means, electronic or mechanical, for any purpose without the express written permission of DDO.